

# Product Sheet

## H\_LILRB2(ILT4) Reporter Jurkat Cell Line


Catalog number: GM-C16968

Version 3.3.1.241206

LILRB2, a member of the leukocyte Ig-like receptor subfamily B (LILRB), is a single-pass transmembrane glycoprotein with extracellular Ig-like domains and intracellular ITIM motifs. Upon binding to ligands, LILRB2 activates ITIM motifs, triggering signals that inhibit T cell activation.

HLA-G, a non-classical MHC class I molecule, is crucial for fetal-maternal tolerance and acts as an immune checkpoint. By interacting with LILRB1 (ILT2) and LILRB2 (ILT4), HLA-G suppresses cytotoxic T cells, NK cells, and B cells, induces T cell anergy, and promotes regulatory T cell (Treg) development. HLA-G on antigen-presenting cells (APCs), such as MDSCs or tolerogenic dendritic cells, further promotes T cell hyporesponsiveness and Treg differentiation.

H\_LILRB2(ILT4) Reporter Jurkat Cell Line is a clonal stable Jurkat cell line constructed using lentiviral technology, constitutive expression of the LILRB2(ILT4) gene, along with signal-dependent expression of a luciferase reporter gene. The reporter cell line is co-cultured with the H\_HLA-G1 OKT3 CHO-K1 Cell Line. The binding of HLA-G1 to LILRB2(ILT4) inhibit T cell signaling. By adding Anti-LILRB2 and Anti-HLA-G1 antibodies, the interaction of HLA-G1 to LILRB2 is blocked, thereby restoring T cell signaling. The luciferase readout indicates the activation level of the signaling pathway, allowing evaluation of the in vitro effects of LILRB1 related drugs.



---

## Specifications

---

<b>Quantity</b>	5E6 Cells per vial, 1 mL
<b>Product Format</b>	1 vial of frozen cells
<b>Shipping</b>	Shipped on dry ice
<b>Storage Conditions</b>	Liquid nitrogen immediately upon receipt

---

<b>Recovery Medium</b>	RPMI 1640+10% FBS+1% P.S
<b>Growth medium</b>	RPMI 1640+10% FBS+1% P.S+0.75 µg/mL Puromycin+3.5 µg/mL Blasticidin
<b>Note</b>	None
<b>Freezing Medium</b>	90% FBS+10% DMSO
<b>Growth properties</b>	Suspension
<b>Growth Conditions</b>	37°C, 5% CO <sub>2</sub>

---

<b>Mycoplasma Testing</b>	The cell line has been screened to confirm the absence of Mycoplasma species.
<b>Safety considerations</b>	Biosafety Level 2
<b>Note</b>	It is recommended to expand the cell culture and store a minimum of 10 vials at an early passage for potential future use.

---

## Materials

---

<b>Reagent</b>	<b>Manufacturer/Catalogue No.</b>
RPMI 1640	VivaCell/C3010-0500
Fetal Bovine Serum	Cegrogen biotech/A0500-3010
Pen/Strep	Thermo/15140-122
Puromycin	Genomeditech/ <a href="#">GM-040401</a>
Blasticidin	Genomeditech/ <a href="#">GM-040404</a>
H_HLA-G1 OKT3 CHO-K1 Cell Line	Genomeditech/ <a href="#">GM-C16834</a>
Anti-H_LILRB2(ILT4) hIgG4 Antibody(MK-4830)	Genomeditech/ <a href="#">GM-27362AB</a>
Anti-H_HLA-G1 hIgG1 Antibody(38373)	Genomeditech/ <a href="#">GM-28208AB</a>
GMOne-Step Luciferase Reporter Gene Assay Kit	Genomeditech/ <a href="#">GM-040503</a>

---

## Figures




Figure 1 | Response to Anti-H\_HLA-G1 hIgG1 Antibody(38373). H\_HLA-G1 OKT3 CHO-K1 Cell Line (Cat. [GM-C16834](#)) was seeded at a density of 2E4 cells per well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H\_HLA-G1 hIgG1 Antibody(38373) (Cat. [GM-28208AB](#)) was incubated with H\_HLA-G1 OKT3 CHO-K1 Cell Line for 1 hour. Subsequently, the H\_LILRB2(ILT4) Reporter Jurkat Cell Line (Cat. GM-C16968) at a concentration of 1E5 cells/well was added, and the coculture proceeded for an additional 5 hours. Firefly luciferase activity is then measured using the GMOne-Step Luciferase Reporter Gene Assay Kit (Cat. [GM-040503](#)). The results indicated maximum blocking folds of approximately [4.5]. Data are shown by drug mass concentration.




Figure 2 | Response to Anti-H\_HLA-G1 hIgG1 Antibody(38373). H\_HLA-G1 OKT3 CHO-K1 Cell Line (Cat. [GM-C16834](#)) was seeded at a density of 2E4 cells per well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H\_HLA-G1 hIgG1 Antibody(38373) (Cat. [GM-28208AB](#)) was incubated with H\_HLA-G1 OKT3 CHO-K1 Cell Line for 1 hour. Subsequently, the H\_LILRB2(ILT4) Reporter Jurkat Cell Line (Cat. GM-C16968) at a concentration of 1E5 cells/well was added, and the coculture proceeded for an additional 5 hours. Firefly luciferase

activity is then measured using the GMOne-Step Luciferase Reporter Gene Assay Kit (Cat. [GM-040503](#)). The results indicated maximum blocking folds of approximately [4.5]. Data are shown by drug molar concentration.




Figure 3 | Response to Anti-H\_LILRB2(ILT4) hIgG4 Antibody(MK-4830). H\_HLA-G1 OKT3 CHO-K1 Cell Line (Cat. [GM-C16834](#)) was seeded at a density of 2E4 cells per well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H\_LILRB2(ILT4) hIgG4 Antibody(MK-4830) (Cat. [GM-27362AB](#)) were incubated with 1E5 cells/well of the H\_LILRB2(ILT4) Reporter Jurkat Cell Line (Cat. [GM-C16968](#)) in a 96-well plate for 1 hour, and then added to the pre-seeded cells. The mixture was incubated for an additional 5 hours. Firefly luciferase activity is then measured using the GMOne-Step Luciferase Reporter Gene Assay Kit (Cat. [GM-040503](#)). The results indicated maximum blocking folds of approximately [11.2]. Data are shown by drug mass concentration.




Figure 4 | Response to Anti-H\_LILRB2(ILT4) hIgG4 Antibody(MK-4830). H\_HLA-G1 OKT3 CHO-K1 Cell Line (Cat. [GM-C16834](#)) was seeded at a density of 2E4 cells per well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H\_LILRB2(ILT4) hIgG4 Antibody(MK-4830) (Cat. [GM-27362AB](#)) were incubated with

1E5 cells/well of the H\_LILRB2(ILT4) Reporter Jurkat Cell Line (Cat. GM-C16968) in a 96-well plate for 1 hour, and then added to the pre-seeded cells. The mixture was incubated for an additional 5 hours. Firefly luciferase activity is then measured using the GMOne-Step Luciferase Reporter Gene Assay Kit (Cat. GM-040503). The results indicated maximum blocking folds of approximately [11.2]. Data are shown by drug molar concentration.




Figure 5 | Response to Anti-H\_LILRB2(ILT4) hIgG4 Antibody(MK-4830). H\_HLA-G1 OKT3 CHO-K1 Cell Line (Cat. GM-C16834) was seeded at a density of 2E4 cells per well in a 96-well plate and incubated overnight. The next day, serial dilutions of the Anti-H\_LILRB2(ILT4) hIgG4 Antibody(MK-4830) (Cat. GM-27362AB) were incubated with 1E5 cells/well of the H\_LILRB2(ILT4) Reporter Jurkat Cell Line (Cat. GM-C16968) CP3 and CP13 in a 96-well plate for 1 hour, and then added to the pre-seeded cells. The mixture was incubated for an additional 5 hours. Firefly luciferase activity is then measured using the GMOne-Step Luciferase Reporter Gene Assay Kit (Cat. GM-040503). Data are shown by drug mass concentration.




Figure 6 | H\_LILRB2(ILT4) Reporter Jurkat Cell Line (Cat. GM-C16968) was determined by flow cytometry using Anti-H\_LILRB2(ILT4) hIgG4 Antibody(1E1) (Cat. GM-27362AB).

## Cell Recovery

Recovery Medium: RPMI 1640+10% FBS+1% P.S

To insure the highest level of viability, thaw the vial and initiate the culture as soon as possible upon receipt. If upon arrival, continued storage of the frozen culture is necessary, it should be stored in liquid nitrogen vapor phase and not at  $-70^{\circ}\text{C}$ . Storage at  $-70^{\circ}\text{C}$  will result in loss of viability.

- a) Thaw the vial by gentle agitation in a  $37^{\circ}\text{C}$  water bath. To reduce the possibility of contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 2 - 3 minutes).
- b) Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by dipping in or spraying with 70% ethanol. All of the operations from this point on should be carried out under strict aseptic conditions.
- c) Transfer the vial contents to a centrifuge tube containing 5.0 mL complete culture medium. And spin at approximately  $176 \times g$  for 5 minutes. Discard supernatant.
- d) Resuspend cell pellet with the recommended complete medium. And dispense the suspension into 1 - 2 T-25 culture flasks.
- e) Incubate the culture at  $37^{\circ}\text{C}$  in a suitable incubator. A 5%  $\text{CO}_2$  in air atmosphere is recommended if using the medium described on this product sheet.

## Cell Freezing

Freezing Medium: 90% FBS+10% DMSO

- a) Centrifuge at  $176 \times g$  for 3 minutes to collect cells.
- b) Resuspend the cells in pre-cooled freezing medium and adjust the cell density to  $5 \times 10^6$  cells/mL.
- c) Aliquot 1 mL into each vial.
- d) Place the vial in a controlled-rate freezing container and store at  $-80^{\circ}\text{C}$  for at least 1 day, then transfer to liquid nitrogen as soon as possible.

## Cell passage

Growth medium: RPMI 1640+10% FBS+1% P.S+0.75  $\mu\text{g}/\text{mL}$  Puromycin+3.5  $\mu\text{g}/\text{mL}$  Blasticidin

Approximately 48-72 hours after the initial thawing, the cells can be passaged for the first time. After this initial passage, the culture medium can be adjusted to growth medium supplemented with antibiotics. If cells are not passaged within 48 hours, it is recommended to add some fresh recovery medium and place the flask horizontally.

- a) When the cell density reaches  $1.5 - 2 \times 10^6$  cells/mL, subculture the cells. Do not allow the cell density to exceed  $2 \times 10^6$  cells/mL.
- b) It is recommended to use T-25 flasks for subculturing.
- c) These cells are suspension cells, and it is recommended to use the "half-medium change" method to maintain optimal cell conditions during passaging.
- d) During passaging, you can directly add fresh growth medium to the culture flask, gently pipette to resuspend the cells, and then transfer the cell suspension to a new T-25 flask for continued culture.



**Subcultivation Ratio: Maintain cultures at a cell concentraion between 3E5 and 1E6 viable cells/mL.**

**Medium Renewal: Every 2 to 3 days**

## Notes

- These cells are sensitive to density, so please ensure that the cell density is maintained within an appropriate range during culture and subculturing.
- During the first passage, pay attention to the nutrient supply; if not subculturing, make sure to add fresh recovery medium every other day as needed.

## Related Products

<b>LILRB2(ILT4)</b>	
<a href="#">H_LILRB2(ILT4) CHO-K1 Cell Line</a>	<a href="#">H_LILRB2(ILT4) HEK-293 Cell Line</a>
<a href="#">Anti-H_LILRB2(ILT4) hIgG4 Antibody(MK-4830)</a>	
<b>LILRB1(ILT2)</b>	
<a href="#">H_LILRB1(ILT2) Reporter Jurkat Cell Line</a>	<a href="#">H_LILRB1(ILT2) CHO-K1 Cell Line</a>
<a href="#">H_LILRB1(ILT2) HEK-293 Cell Line</a>	<a href="#">Rhesus_LILRB1 CHO-K1 Cell Line</a>
<a href="#">Anti-LILRB1(ILT2) mIgG1 Antibody(12D12)</a>	
<b>LILRB4(ILT3)</b>	
<a href="#">H_LILRB4(ILT3) CHO-K1 Cell Line</a>	<a href="#">H_LILRB4(ILT3) HEK-293 Cell Line</a>
<b>LILRB5</b>	
<a href="#">H_LILRB5 CHO-K1 Cell Line</a>	
<b>LILRB3(ILT5)</b>	
<a href="#">H_LILRB3 Reporter Jurkat Cell Line</a>	<a href="#">H_LILRB3(ILT5) CHO-K1 Cell Line</a>
<a href="#">H_LILRB3(ILT5) HEK-293 Cell Line</a>	
<a href="#">Anti-LILRB3 hIgG1 Antibody(7C5)</a>	

## Limited Use License Agreement

Genomeditech (Shanghai) Co., Ltd grants to the Licensee all intellectual property rights, exclusive, non-transferable, and non-sublicensable rights of the Licensed Materials; Genomeditech (Shanghai) Co., Ltd will retain ownership of the Licensed Materials, cell line history packages, progeny, and the Licensed Materials including modified materials.

Between Genomeditech (Shanghai) Co., Ltd, and Licensee, Licensee is not permitted to modify cell lines in any way. The Licensee shall not share, distribute, sell, sublicense, or otherwise provide the Licensed Materials, or progenitors to third parties such as laboratories, departments, research institutions, hospitals, universities, or biotechnology companies for use other than for the purpose of outsourcing the Licensee's research.

Please refer to the Genomeditech Cell Line License Agreement for details.